


12.141

Electron Microprobe Analysis



Nilanjan Chatterjee, Ph.D.

Principal Research Scientist

Electron Microprobe Facility

Department of Earth, Atmospheric and Planetary Sciences



Electron Microprobe or Electron Probe Microanalyzer (EPMA)

Surface characterization of solids at the micrometer-scale:

- Surface topography and compositional imaging
- Complete chemical analysis of microscopic volumes
(Electron beam induced X-ray emission spectrometry)


Signals produced in the Electron Microprobe

*Electron
beam*

Cathodoluminescence (CL)

Back-scattered electron (BSE)

Characteristic X-ray



Secondary
electron (SE)

Specimen

Qualitative analysis

- Visual characterization and identification of phases in image (shape, size, surface relief, etc.)
- Identification of elements in each phase (no concentration measurement)

Semi-quantitative analysis

- Spatial distribution of elements in an image
- Quick and approximate spot concentration measurement without calibration

Quantitative analysis


- Full quantitative micro-chemical analysis

Concentration of all elements present at the spot

- Elemental concentration mapping


Concentration of all elements present at each pixel of the image

Electron-specimen interactions



Elastic Scattering


- Back-scattered electron




Inelastic Scattering

- Characteristic X-rays
- Secondary electron
- Cathodoluminescence

Electron interaction volume



Electron interaction volume




$$R = 0.0276 E^{1.67} \frac{A}{\rho Z^{0.889}}$$

(A = atomic weight, ρ = density)

(Kanaya-Okayama range)

- Increases with electron beam energy, E
- Decreases with sample atomic number, Z

Electron interaction volume





Typical ranges (15 kV, perpendicular beam):

C ($Z = 6$) **1.8 μm**

Fe ($Z = 26$) **1.1 μm**


U ($Z = 92$) **0.8 μm**

Electron probe diameter and Electron interaction volume



Elastic scattering cross-section

For a scattering angle $> \phi_e$, cross-section
(events.cm²/e⁻.atom)



$E_1 = E_0$, large ϕ_e

$$Q_e = 1.62 \times 10^{-20} (Z^2/E^2) \cot^2(\phi_e/2)$$

Z : atomic number


E : beam energy

ϕ_e : elastic scattering angle

- Increases with sample atomic number, Z
- Decreases with electron beam energy, E


Back-scattered electron (BSE)

(Elastic scattering)



- Beam electrons scattered backward from specimen surface
- High energy electrons with energy about the same as that of the electron beam
- **BSE image resolution improves with shrinking of the electron interaction volume through:**
 1. Decrease in beam energy
 2. Increase in specimen atomic number


Electron backscatter coefficient



BSE image contrast is better among low Atomic Number elements

Backscattered electron image


Back-scattered electron



Polished surface

Function of
composition


Plane polarized transmitted light







Thin section

Function of optical
properties

The X-ray spectrum





Phase identification: EDS X-ray spectra




Z	1	8	11	12	13	14	19	20	22	26
ilm:		O							Ti	Fe
grt:		O		Mg	Al	Si		Ca		Fe
bt:	H	O		Mg	Al	Si	K			Fe
hbl:	H	O		Mg	Al	Si		Ca		Fe
plg:		O	Na		Al	Si		Ca		

JEOL JXA-8200 Superprobe




Sample surface is perpendicular to the electron beam

Back-scattered electron detector




Located vertically above the specimen
A split ring shape



solid-state diode



Compositional and topographic imaging with BSE detector



A+B
Compositional
mode

A-B
Topographic
mode


Energy Dispersive X-ray Spectrometer (EDS)



- *EDS detector: a 'p-n' layer of intrinsic **Si(Li)** semiconductor; Be window; aperture wheel*
- *Multichannel analyzer (MCA) processes the X-ray signal*

Secondary electron (SE)

(Inelastic scattering)





$E_1 < E_0$, small ϕ_i


- Electrons from specimen surface are mobilized by beam electrons
- Emitted at low energies (typical: <10 eV)

(recall BSE are high energy beam electrons that underwent elastic scattering, $E_1 = E_0$, E_0 typically being 10-20 keV)


Secondary electron detector



Located on the side wall of the sample chamber




Imaging with the E-T detector



-ve Faraday cage bias

only BSE

Surfaces in direct line of sight
are illuminated




+ve Faraday cage bias

BSE + SE

All surfaces are illuminated

Cathodoluminescence (CL)

Light generated from semiconductor samples through electron beam interaction



- *Pure material has an empty conduction band; does not conduct*
- *Trace impurities add additional energy levels in the band gap that can accept electrons in the excited state*
- *CL photon is emitted as electron drops back to the valence band*

Cathodoluminescence spectrometer

Optical microscope camera (not used)





Optical microscope light (turned off)



Optical spectrometer

Cathodoluminescence spectrometry



MIT OpenCourseWare
<http://ocw.mit.edu>

12.141 Electron Microprobe Analysis
January (IAP) 2012

For information about citing these materials or our Terms of Use, visit: <http://ocw.mit.edu/terms>.